



In situ X-ray for monitoring crystallization and phase transitions in suspension: In-situX®

Simon Clevers, A. Lafontaine, G. Baaklini, M. Schindler, Morgane Sanselme,
Yohann Cartigny, V. Dupray, Gérard Coquerel

► To cite this version:

Simon Clevers, A. Lafontaine, G. Baaklini, M. Schindler, Morgane Sanselme, et al.. In situ X-ray for monitoring crystallization and phase transitions in suspension: In-situX®. BIWIC 2019 International Workshop on Industrial Crystallization, Aug 2019, Rayong, Thailand. . hal-02295903

HAL Id: hal-02295903

<https://normandie-univ.hal.science/hal-02295903>

Submitted on 24 Sep 2019

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



In situ X-ray for monitoring crystallization and phase transitions in suspension: In-situX®




S. Clevers, A. Lafontaine, G. Baaklini, M. Schindler, M. Sanselme, Y. Cartigny, V. Dupray and G. Coquerel

Laboratoire Sciences et Méthodes Séparatives, EA3233, Université de Rouen Normandie,
76821 Mont Saint Aignan, France.
valerie.dupray@univ-rouen.fr



Introduction




Not suitable for:

- In-situ analyses
- Metastable phase
- Efflorescent solvate
- Crystals in suspension
- Transient states

Aim of the study

To highlight, through different examples, the possibilities of In-situX® as analytical tool able to solve academic as well industrial problems

Conclusion




- No need for sampling (cytotoxic compounds, radioactive elements, low temperature crystallization, efflorescent solvates, CO_2 , O_2 , light, moisture sensitive compounds)
- Laboratory X-ray Source

Solid – solid transition studies of all sorts:


- Polymorphism (kinetics of transition during slurring)
- De/resolvation (peritectic transitions in binary, ternary systems, etc.)
- Co-crystal formation
- Host Guest associations
- Clathrates

- Optimization of the resolution of chiral molecules
- Optimization of the filterability by means of temperature cycling


DMU/H₂O- Polymorphic transition




Detection of solvate




Determination of peritectic invariants → Non-solvated phase can be obtained directly in suspension



Vapor/solid equilibria



In-situX® - Principle




Diffractometer[5]:

- Geometry: θ - θ
- Detector: Lynx eyes™
- X-ray tube: anticathode in copper

Computer-controlled (stirring rate, chemical injection, temperature) by Labview® software

Measure in suspension

⇒ $V = 40\text{mL}$, $V_{\max} = 100\text{mL}$




Vapor/solid equilibria

⇒ Measurement vs. relative humidity

AFRODISM[6]:

- Temperature control of the exterior of the reactor
- Anti-frosting and anti-condensation on the Kapton window



References

- [1] C. Näther, C. Döring, I. Jess, P. G. Jones, and C. Taouss, *Acta Crystallogr B Struct Sci*, vol. 69, no. 1, 2012.
- [2] G. Baaklini *et al.* "Depression of the polymorphic transition temperature of 1,3-dimethylurea (DMU) by insertion of hundreds of ppm of water in the crystal lattices and crystal structure determination of DMU monohydrate." to be published
- [3] A. Lafontaine, PhD Thesis, Université de Rouen , 2013
- [4] A. Lafontaine, M. Sanselme, Y. Cartigny, P. Cardinael, and G. Coquerel, *Journal of Thermal Analysis*, vol. 112, no. 1, pp. 307–315, Nov. 2012.
- [5] G. Coquerel, M. Sanselme, and A. Lafontaine, US9255898B2 patent, 2014
- [6] G. Coquerel, S. Clevers EP18305580.5 patent, 2018